Indication and Important Safety Information

gammaCore Sapphire™ (non-invasive vagus nerve stimulator) is intended to provide non-invasive vagus nerve stimulation (nVNS) on the side of the neck. gammaCore is indicated for:

- Preventive treatment of cluster headache in adults
- Acute treatment of cluster headache in adults
- Preventive treatment of migraine headache in adolescents (age 12 and older) and adults
- Acute treatment of migraine headache in adolescents (age 12 and older) and adults

The long-term effects of the chronic use of the device have not been evaluated.

Safety and efficacy of gammaCore has not been evaluated and therefore is NOT indicated for:

- Adolescents with congenital cardiac issues

gammaCore contraindications include but are not limited to:

- Patients with an active implantable medical device, such as a pacemaker, hearing aid implant, or any implanted electronic device
- Patients who are using another device at the same time (e.g., TENS Unit, muscle stimulator) or any portable electronic device (e.g., mobile phone)

Warnings

- Safety and efficacy of gammaCore have not been evaluated in the following patients:
 - Patients diagnosed with narrowing of the arteries (carotid atherosclerosis)
 - Patients with a metallic device such as a stent, bone plate or bone screw implanted at or near their neck
 - Patients who have had surgery to cut the vagus nerve in the neck (cervical vagotomy)
 - Pediatric patients (less than 12 years of age)
 - Pregnant women
 - Patients with clinically significant hypertension, hypotension, bradycardia, or tachycardia

Precautions

- The long-term safety and effectiveness of the gammaCore device has not been demonstrated in adolescents 12-17 years of age. Due to hormonal and cognitive development changes in adolescents, this population should be closely monitored while using the device. The use of the device in this population is based on extrapolated data from a clinical study in adults.
- You must read the gammaCore *Instructions for Use* before using gammaCore
- Only use gammaCore as described in these *Instructions for Use* or as otherwise directed by your Healthcare Provider
- Only use an electroCore-approved electrode gel with gammaCore.

Note: This list is not all inclusive. Please refer to the gammaCore Instructions for Use for all of the important contraindications, warnings and precautions before using or prescribing this product.

Healthcare Provider authorization required.

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