## Warnings and Precautions - Important Safety Information gammaCore™ (non-invasive vagus nerve stimulator)

**Warnings indicate instructions, which, if not followed, may result in serious injury or death to the device user or to the patient.**

**Precautions indicate instructions, which, if not followed, may result in damage to the equipment or degradation in the quality of treatment.**

### Indications

gammaCore™ is a medical device, non-invasive vagus nerve stimulator (nVNS) that is intended to provide mild electrical stimulation to the vagus nerve on the side of the neck.

gammaCore is indicated for:

- Preventive treatment of cluster headache in adults
- Acute treatment of cluster headache in adults
- Preventive treatment of migraine headache in adolescents (age 12 and older) and adults
- Acute treatment of migraine headache in adolescents (age 12 and older) and adults

The main benefits associated with the use of gammaCore are a reduction in headache days per month or attack frequency and reduced pain.

### Warnings

- The long-term effects of the chronic use of gammaCore have not been evaluated.
- Safety and efficacy of gammaCore has not been evaluated and therefore is NOT indicated for:
  - Adolescents with congenital cardiac issues
- Safety and efficacy of gammaCore have not been evaluated in the following patients:
  - Patients with uncontrolled hypertension, hypotension, bradycardia, or tachycardia
  - Patients with a history of baseline cardiac disease or atherosclerotic cardiovascular disease, including congestive heart failure, known severe coronary artery disease, or recent myocardial infarction (within 5 years)
  - Patients with a metallic device such as a stent, bone plate or bone screw implanted at or near their neck
  - Patients with a history of abnormal baseline ECG, prolonged QT interval or arrhythmia
  - Patients who have had surgery to cut the vagus nerve in the neck (cervical vagotomy)
  - Pediatric patients
  - Pregnant women
  - Patients with active cancer or cancer in remission
  - Patients with an abnormal cervical anatomy
  - Patients with a history of brain tumor, aneurysm, bleed, or head trauma
  - Patients with a history of syncope or seizures
- Contact your HealthCare Provider (HCP) if your symptoms continue or worsen
- Treatment is intended to be given (administered) as directed by an HCP. Your HCP or Customer Service must train you in the proper use of gammaCore

### Precautions

- The long-term safety and effectiveness of gammaCore have not been demonstrated in adolescents 12-17 years of age. Due to hormonal and cognitive development changes in adolescents, this population should be closely monitored while using the device. The use of the device in this population is based on extrapolated data from a

clinical study in adults.

***Before Use:***

- Only use gammaCore as described in the User Manual, or as otherwise directed by the HCP
- The use of more than 24 stimulations per day has not been evaluated in controlled clinical trials. Do not use more than 24 stimulations in a 24-hour period
- Only use an electroCore-approved gel with gammaCore. Please contact Customer Service for an electroCore-approved gel that works with the device
- Do not share gammaCore with another person
- gammaCore should not be applied across or through the head, directly on the eyes, covering the mouth, on the chest or the upper back, or over the heart
- Remove jewelry that may touch the treatment location (necklaces, earrings, etc.)
- Carefully examine the device for any signs of damage or defects

***Do not use gammaCore if:***

- The stimulation surfaces are broken or cracked
- The casing is cracked, dented, or appears to be damaged
- An “E” with a numeric code is displayed on the screen when the device is turned on. This means that there is an error

***During Treatment:***

***Discontinue treatment if you experience:***

- Light-headedness, dizziness, or chest pain
- Excessive skin irritation

If the device seems to malfunction, discontinue use immediately, continue taking usual medications and seek medical care as needed. When possible, contact Customer Service for assistance with your device; Customer Service cannot provide medical assistance.

***Caring for Your Device:***

- Keep gammaCore away from water or other liquids, including cleaning liquids.
- Moisture may damage the device. Keep gammaCore away from items such as nebulizers and steam kettles.
- Store gammaCore in a safe location out of reach of children.
- Avoid exposure to extreme hot or cold temperatures outside the range of 0°C to 38°C (32°F to 100°F). Exposure to such conditions may cause the device to malfunction.
- Do not attempt to replace the device battery. If the device is not working, contact Customer Service.
- Do not open or take apart the case or attempt to repair or modify the device. There are no user serviceable parts. If the device is not working, contact Customer Service.
- As the device contains a lithium-ion battery, do not intentionally damage, burn, or puncture the device.
- Wireless communications equipment, such as wireless home network devices, mobile phones, cordless telephones and their base stations, and walkie-talkies, can affect this equipment. Keep gammaCore at least 3.3 m (10.8 ft) away from these items while in use.

**Contraindications**

Contraindications include:

- Patients with an active implantable medical device, such as a pacemaker, hearing aid implant, or any implanted electronic device
- Patients who are using another device at the same time (e.g., TENS Unit, muscle stimulator) or any portable electronic device (e.g., mobile phone)

Patients should not use gammaCore:

- While driving, operating machinery, or during any activity that may put patient at risk of injury
- Near microwave machines, magnetic resonance imaging, radio frequency surgical, or computer-aided tomography machines
- In an explosive atmosphere or in the presence of flammable gas mixtures
- If patient has an open wound, rash, infection, swelling, cut, sore, drug patch, or surgical scar(s) on the neck at the treatment location
- If patient has wet skin, is in the water, or just stepped out of the water (e.g., shower, bath, pool)

## **Safety and Performance**

Please visit [www.gammacore.com/prescribing-gammacore/clinical-efficacy](http://www.gammacore.com/prescribing-gammacore/clinical-efficacy) for additional Safety and Performance information.

**Caution: Healthcare Provider (HCP) Authorization Required.**

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